

food tests

FOR ILLITERATE MONKEYS



AWARDED THE
BEST COOK
BOOK ON BBC

written by

kesavan karadiar




published by

nasim haque





what to GET



Spatula 1

Pipette 1

Dish 1

Measuring cylinder 1


gloves 1 pair

safety goggles 1 pair

lab coat 1

iodine solution >5 ml

food sample >1 cm³





what to DO

- 1 Place a small amount of solid food in a dish, or pour 1 cm³ of liquid food into a measuring cylinder.**
- 2 Use a pipette to draw some iodine solution.**
- 3 Add 3 to 4 drops of iodine onto the food sample.**
- 4 Wait a few seconds and observe the colour change.**
- 5 If the solution turns blue or black, starch is present; if not, starch is absent.**



what to GET



test tube 1

test tube rack 1

measuring cylinder 1

Benedict's solution 5 ml

pipette 1



food sample $>1\text{cm}^3$

spatula 1

goggles 1 pair

lab coat 1

gloves 1 pair





what to DO

- 1 Place a small amount of the food sample in a test tube. If using a liquid food, pour about 1 cm³ into the tube.**
- 2 Add an equal volume of Benedict's solution to the food sample.**
- 3 Gently heat the test tube in a hot water bath for about 2-3 minutes.**
- 4 Remove the tube and observe the colour change.**
- 5 If the solution turns green, yellow, orange, or brick-red, reducing sugars are present. If it stays blue, no reducing sugars are present.**



what to GET



test tube 1

test tube rack 1

measuring cylinder 1

Biuret solution 5 ml

pipette 1



food sample $>1\text{cm}^3$

spatula 1

goggles 1 pair

lab coat 1

gloves 1 pair





what to DO

- 1 Place a small amount of the food sample in a test tube. If using a liquid food, pour about 1 cm³ into the tube.**
- 2 Add an equal volume of Biuret solution to the food sample.**
- 3 Gently shake the test tube to mix the contents.**
- 4 Wait a few seconds and observe the colour change.**
- 5 If the solution turns purple or violet, protein is present; if it stays blue, protein is absent.**

IS IT CAKE?



what to GET



test tube 1

test tube rack 1

measuring cylinder 1

ethanol 5 ml

pipette 1



food sample $>1\text{cm}^3$

spatula 1

goggles 1 pair

lab coat 1

gloves 1 pair



what to DO

- 1 Place a small amount of the food sample in a test tube. If using a solid, crush it first.**
- 2 Add about 2 cm³ of ethanol to the test tube.**
- 3 Shake the tube gently for about 30 seconds to dissolve any lipids.**
- 4 Add 2 cm³ of cold distilled water to the mixture.**
- 5 Observe the solution. If a cloudy white emulsion forms, lipids are present. If the solution remains clear, lipids are absent.**